

## OXYGEN ADMINISTRATION

Start Time: \_\_\_\_\_

Stop Time: \_\_\_\_\_ Date: \_\_\_\_\_

Candidate's Name: \_\_\_\_\_

Evaluator's Name: \_\_\_\_\_

	Points Possible	Points Awarded
Takes, or verbalizes, body substance isolation precautions	1	
Assembles the regulator to the tank	1	
Opens the tank	1	
Checks for leaks	1	
Checks tank pressure	1	
Attaches non-rebreather mask to oxygen	1	
Prefills reservoir	1	
Adjusts liter flow to 12 liters per minute or greater	1	
Applies and adjusts the mask to the patient's face	1	
<b>Note: The examiner must advise the candidate that the patient is not tolerating the non-rebreather mask. The medical director has ordered you to apply a nasal cannula to the patient.</b>		
Attaches nasal cannula to oxygen	1	
Adjusts liter flow to six (6) liters per minute or less	1	
Applies nasal cannula to the patient	1	
<b>Note: The examiner must advise the candidate to discontinue oxygen therapy</b>		
Removes the nasal cannula from the patient	1	
Shuts off the regulator	1	
Relieves the pressure within the regulator	1	
<b>Total :</b>	<b>15</b>	

### Critical Criteria

\_\_\_\_\_ Did not take, or verbalize, body substance isolation precautions.

\_\_\_\_\_ Did not assemble the tank and regulator without leaks.

\_\_\_\_\_ Did not prefill the reservoir bag.

\_\_\_\_\_ Did not adjust the device to the correct liter flow for the non-rebreather mask.  
(12 liters per minute or greater)

\_\_\_\_\_ Did not adjust the device to the correct liter flow for the nasal cannula.  
(6 liters per minute or less)